

# 2024 Clinical Recommendations for the Treatment of Unilateral Hearing Loss/Single-Sided Deafness with Cochlear Implantation

In collaboration with participants of the Cochlear sponsored UHL/SSD Thought Leadership Conferences

## Introduction

Single-sided deafness (SSD) impacts over 345,000 adults in the United States, with an estimated 5.2% experiencing some degree of unilateral hearing loss (Kay-Rivest et al, 2021). SSD significantly impacts spatial hearing, speech understanding in noise, and quality of life, affecting various auditory and non-auditory functions (Dillon et al., 2017b; Kumpik & King, 2019). Historically, treatment options for SSD have been limited. Many patients have gone untreated or relied on treatment options which stimulate the contralateral ear, such as contralateral routing of signal (CROS) hearing aids or bone conduction solutions. Today, expanding criteria in conjunction with advancements in hearing technologies, such as cochlear implants (CIs) and bone conduction devices, have significantly changed the treatment landscape for individuals with SSD. Since approval in early 2022, SSD accounts for between 13-16% of CI registrations. This number continues to grow, as cochlear implantation is the only solution which allows for binaural input from two ears, providing potential benefits and improved outcomes over other treatment options.

Early experiences with SSD have sparked new investigations into individual patient factors, candidacy criteria, and incidence of non-use. These insights have shifted the focus during implant candidate selection with increased emphasis on quality of life, expectations, and motivation.

This paper will provide comprehensive information on the treatment of SSD with a CI, based on recommendations from industry experts and Cochlear North America, two years post-FDA approval. Key areas of focus include:

- Adult treatment pathway and candidacy considerations
- Patient selection criteria, quality of life assessment, patient motivation, and establishing realistic expectations for performance and rehabilitation
- Essential aspects of programming a cochlear implant for patients with SSD
- Follow-up protocols to ensure optimal outcomes for adults

A list of professionals supporting this consensus can be found in Appendix A.

## SSD thought leadership workshops


In preparation for FDA approval of CI for SSD, an inaugural Thought Leadership Workshop was held in December 2021. This collaborative effort led to the development of a comprehensive SSD treatment pathway and development of a white paper focused on quality-of-life impacts, medical history, audiological profile, imaging, device programming, and post-implant assessments for adults. A follow-up Thought Leadership Workshop was conducted in May 2024. Moderated by Doctors Kevin Brown and Daniel Zeitler, the agenda included a mix of scientific presentations and panel discussions, aimed at refining and enhancing candidacy and treatment protocols based on the latest clinical experiences and outcomes. The critical updates to the previous guidance for adults include the following:

- Counseling during the candidacy evaluation phase is critical and assessment of patient expectations, motivation, and hearing deficits should be standard practice. Additionally, it is crucial to establish the expectation that outcomes are contingent upon device usage and aural rehabilitation.
- The SSD pre- and post-implant evaluation protocol was updated to align with the Minimum Speech Test Battery, Version 3 (MSTB-3), utilizing AzBio Sentences with a fixed signal-to-noise ratio (SNR) as the primary metric for evaluating speech understanding in noise, as an alternative to the Bamford-Kowal-Bench Speech-in-Noise (BKB-SIN).
- The post-implant activation and follow-up schedule should align with current care model practices and be updated to include intervals at Initial Activation, 1 month, 3 months, and 12 months. This eliminates the 6-month visit from the previous treatment pathway.

## Navigating single-sided deafness: Adult treatment pathway

The Adult Treatment Pathway (Figure 1.) summarizes the steps in the evaluation process including important considerations and assessments which help the implant team determine the most appropriate intervention and guides counseling. Refer to Appendices B and C for additional references, resources, and supporting discussion.

Figure 1: Adult Treatment Pathway for Single-Sided Deafness



## Indications

The Cochlear™ Nucleus® Implant system indication was expanded in January 2022 to include patients with SSD aged 5 years or older. However, depending on patient history, etiology of hearing loss, and other factors, bone conduction solutions, non-surgical options, or no treatment/observation may be considered in some cases. The current Nucleus® Cochlear Implant and bone conduction Baha® and Osia® systems indications are listed below in Figures 2 and 3. There is overlap in the indications for these implantable solutions which provides the clinician with a range of treatment options available to meet individual patient needs.

Figure 2.

### Cochlear Implants for UHL/SSD 5 years of age or older

| Ear to be Implanted   | Contralateral Ear   |
|---|---|
| Severe to profound sensorineural hearing loss defined as:<br>Pure-tone average at .5, 1, 2, 4 KHz >80 dB HL | Normal or near normal hearing defined as:<br>Pure-tone average at .5, 1, 2, 4 kHz <30 dB HL                                 |
| Aided CNC word score or developmentally appropriate<br>word test <5%  | Recommended 2 weeks to 1 month experience wearing<br>appropriately fit CROS hearing aid or other<br>suitable hearing device |

Contraindications for cochlear implantation include: 1) absence of cochlea development, 2) absence of cochlear nerve, 3) active middle ear infections, 4) tympanic membrane perforation in the presence of active middle ear disease, or 5) duration of profound sensorineural hearing loss greater than ten years.

Figure 3.

### Bone Conduction for UHL/SSD 5 years of age or older

| Poorer Hearing Ear                  | Contralateral Ear   |
|-------------------------------------|---|
| Profound sensorineural hearing loss | Normal or near normal hearing defined as: Pure-tone average at<br>.5, 1, 2, 3 kHz <= 20 dB HL |

Contraindications for bone conduction implantation include the following: 1) insufficient bone quality or quantity to support implantation of both the BI300 Implant and the OSI300 Implant, 2) chronic or non-revisable vestibular or balance disorders that could prevent benefit from the device, as determined by good clinical judgment, 3) abnormally progressive hearing loss or evidence that hearing loss is bilateral, retrocochlear or bilateral central origin, 4) evidence of conditions that would prevent good speech recognition potential as determined by good clinical judgment, or 5) skin or scalp conditions that may preclude attachment of the Sound Processor or that may interfere with the use of the Sound Processor.

---

## Candidacy evaluation

This section provides the professional with a suggested protocol to evaluate and assess hearing performance in the SSD patient.

### Medical Evaluation

For individuals with SSD, a typical medical evaluation includes an extensive otologic history with careful attention to details regarding chronic ear disease, history of vestibular dysfunction, previous ear or skull-based surgery, meningitis, bacterial labyrinthitis, autoimmune inner ear disease other comorbidities. A thorough physical examination including a cranial nerve exam and microscopic otoscopy should also be performed. Etiology of the hearing loss with consideration for risk of progression in the affected ear or contralateral ear is also important to consider. In addition to this, tinnitus is prevalent in this population and can significantly impact quality of life, making it crucial to incorporate tinnitus into the overall treatment plan and evaluation process. Imaging is recommended due to the asymmetry of the hearing loss and to confirm cochlear structural integrity in cases of congenital hearing loss. An MRI will document the anatomy and patency of the cochlea as well as allow visualization of the internal auditory canal to provide insight into possible cochlear nerve deficiencies or retrocochlear pathology which could impact potential outcomes. Duration and progression of the hearing loss should also be considered. Shorter duration of deafness (severe-to-profound hearing loss < 10 years) may result in improved outcomes as suggested by Rahne & Plontke (2016). However, other studies suggest duration of deafness for SSD has less impact on outcomes compared to duration in cases of bilateral sensorineural hearing loss (Nassiri et al. 2021a). Thus duration alone should not preclude consideration and should be considered along with other patient variables.

### Audiological Evaluation

The pre-operative candidacy evaluation should establish baseline measurements for both objective and subjective measures to confirm hearing loss meets FDA or institutional audiometric criteria, demonstrate the need for intervention and provide a basis for pre- and post-implant comparison. This includes routine audiometric testing and evaluation of the speech understanding abilities with the better ear either plugged and/or masked depending on patient age or risk of central inhibition.

Specifically, evaluation of aided and unaided CNC words, as well as AzBio Sentences in noise in the ear to be implanted (65 dBA with 60-65 dB of masking to the normal hearing ear) is recommended to align with MSTB-3 recommendations. Aided binaural testing is not required for SSD candidates; however, confirmation of normal hearing levels in the contralateral ear is necessary for reimbursement purposes.

For additional information regarding MSTB-3 materials and set-up refer to <https://www.cochlearimplanttraining.com/minimum-speech-test-battery>.

### Assessment of Quality of Life, Goals, Motivation and Expectations

Emphasis should be placed on using one or more validated objective assessments to better understand individual impact on quality of life, tinnitus handicap, listening fatigue, and other reported challenges. Several tools are available to evaluate and understand an individual's unique needs. The Cochlear Implant Quality of Life Profile (CIQOL), developed by McRacken and Hand (2019), offers insights into a patient's functional hearing ability. The Speech, Spatial, and Qualities of Hearing Scale (SSQ-12), created by Noble et al. (2019), provides valuable information about spatial hearing in real-world scenarios, particularly beneficial for the SSD population. The Tinnitus Handicap Inventory (THI) by Newman et al. (1996) should be considered when tinnitus is reported. All these measures can be administered both pre- and post-implantation to monitor progress. These tools can also be used during the counseling process to evaluate and discuss patient and family expectations and establish goals. This is critical, as alignment with realistic outcome expectations is often tied to patient satisfaction. Appendix C provides additional description of these assessments, the information they offer, and how they can be administered.

Additionally, thought leaders agree evaluating patient and family motivation to pursue treatment is essential, along with their understanding and willingness to complete aural rehabilitation. Expectations for device wear time of 10+ hours daily and regular direct streaming should also be established.

## Treatment options

When discussing treatment options with SSD candidates, the importance of understanding quality of life impacts, as well as establishing expectations and goals has become evident. Understanding patient motivation for treatment and other co-existing complaints, such as tinnitus, is also critical and may help guide discussion. Observation might be considered for patients who are not interested in treatment or do not perceive a quality-of-life impact as a result of their hearing loss. Non-surgical options such as CROS or Baha SoundArc™ have historically been considered in cases where the patient is seeking treatment but medical comorbidities or other factors preclude surgical intervention. Bone conduction implants such as Osia systems offer viable solutions in cases where structural or neural integrity is a concern or when the hearing history does not support CI candidacy. Patients using Baha and Osia systems report numerous benefits, including improved speech in noise perception and improvements in quality of life (Almugathwi et al., 2020). However, both non-surgical and bone conduction solutions rely on transcranial re-routing of sound, which limits more complex binaural sound processing capability which can only be achieved with a cochlear implant. Presence of tinnitus may also influence treatment options since cochlear implants have been shown to significantly reduce tinnitus in a substantial 50-90% of cases with only a small percentage experiencing no change or worsening of their tinnitus (Idriss et al, 2022, Assouly et al, 2021, Borges et al, 2021, and Levy et al, 2019). Treatment options and potential auditory benefits have been summarized below (Figure 4).

Table 4 : Potential Auditory Benefit of Treatment Options

|                               |                         | TREATMENT OPTION             |             |                       |                              |
|-------------------------------|-------------------------|------------------------------|-------------|-----------------------|------------------------------|
|                               |                         | No Treatment/<br>Observation | CROS/BiCROS | Baha®/Osia®<br>System | Nucleus®<br>Cochlear Implant |
| POTENTIAL<br>AUDITORY BENEFIT | Overcome Head<br>Shadow | ×                            | ✓           | ✓                     | ✓                            |
|                               | Sound Lateralization    | ×                            | ✓           | ✓                     | ✓                            |
|                               | Improved Localization   | ×                            | ×           | ×                     | ✓                            |
|                               | Binaural Summation      | ×                            | ×           | ×                     | ✓                            |
|                               | Squelch                 | ×                            | ×           | ×                     | ✓                            |

---

## Post-Implant programming & assessment

### Device Programming for SSD/UHL

Due to the presence of a normal hearing ear, the following recommendations are for all SSD/UHL programming sessions:

- Sound from the computer should be disabled/ muted during the programming session
- Use of an earplug to minimize any cues of stimulation via the normal hearing ear
- Minimize background noise in the test environment
- Limit conversation during programming

Population mean is the standard starting point at initial activation and for measurement of psychophysics. The use of population mean enables the audiologist to easily create a MAP that is audible and comfortable without needing to make numerous threshold (T) and comfort (C) level measurements during the initial activation. This method creates a starting MAP using mean T- and C-level profiles derived from a large global patient data set to generate a starting MAP (Cochlear Limited 2020) which is electrode specific. A population mean MAP starts with a dynamic range of 46 clinical levels (CLs) and typically has thresholds of 70-80 CLs and comfort levels of 110-120 CLs. The T- and C-level profile is gradually increased while maintaining the dynamic range until the patient perceives the sound to be “loud.” Sweeping using banded C-levels can then be performed to account for any profile differences.

Since the SSD population has a natural comparison between normal hearing and hearing with a CI, there may be some sound quality adjustments that need to be addressed. Mismatch between ears may be resolved by making changes to the frequency allocation table (FAT) within Custom Sound® Pro. To optimize sound quality in these cases, MAPs with the default and alternative frequency allocation table assignments are recommended (Landsberger et al, 2015, Sagi et al, 2021, and Tóth et al, 2023). MAPs should be created with a low frequency cut-off of 188 Hz (default), 438 Hz and 563 Hz and trialed at initial activation. Patients should be encouraged to try each of the MAPs to assess sound quality prior to returning for their first follow-up visit.

Clinical experience collected following FDA approval supports this approach as a simple and effective way to overcome subjective mismatch between ears and does not require additional imaging.

At follow-up programming sessions, further optimization of the MAP and confirmation of audibility should occur. T-levels should be measured at least one time. With experience, recipients should be able to provide consistent measurement of sound field detection levels to confirm T-levels are set appropriately and to confirm good audibility of sound (detection between 20-30 dB HL) with consideration to the frequency allocation utilized in the MAP. The low frequency boundary would be expected to impact detection at 250 and/or 500 Hz, depending on the settings. Confirmation of these levels should be conducted via sound-field audiograms in the booth with the better ear plugged, or via direct streaming using Remote Check. T-levels would only be remeasured if indicated by sound-field testing or remote check audiograms. C-levels can also be optimized and adjusted based upon Electrically Evoked Stapedial Reflex Thresholds (ESRT). Neural Response Telemetry (NRT) is another tool which may be utilized to help confirm mapping.

The following evaluation protocol, based on Thought Leadership collaboration, outlines recommended assessments and intervals (Figure 5).



Figure 5.

| Initial Activation  | 1-Month Visit<br>optional   | 3 and 12-Month Visits   |
|---|---|---|
| New MAP created based on Population Mean  | Determine preferred frequency allocation  | Aided soundfield thresholds   |
| Provide MAPs with alternative frequency allocations (188, 438, and 563 Hz)                                | Aided soundfield thresholds   | <b>AzBio Sentences in noise</b><br>(65 dBA with 60-65 dB masking to normal ear)     |
| <b>Counseling</b><br>Datalogging 10+ hours of wear time/day<br>Review preoperative goals and expectations | <b>AzBio Sentences in noise</b><br>(65 dBA with 60-65 dB masking to normal ear)     | <b>Optimize T/C-levels</b><br>(As needed based on detection and performance scores) |
| <b>Rehabilitation</b><br>Direct streaming for 30-60 minutes per day, 5 days/week                          | <b>Optimize T/C-levels</b><br>(As needed based on detection and performance scores) | <b>Counseling</b><br>Datalogging 10+ hours of wear time/day                         |
|   | <b>Counseling</b><br>Datalogging 10+ hours of wear time/day                         | <b>Rehabilitation</b><br>Direct streaming for 30-60 minutes per day, 5 days/week    |
|   | <b>Rehabilitation</b><br>Direct streaming for 30-60 minutes per day, 5 days/week    | Questionnaire(s)<br>(Compare to preop)  |

**Remote Check allows for complete isolation of the implant ear for detection thresholds, hearing performance in noise and monitoring of datalogging. Remote Check should be considered to augment the in-person assessment or in lieu of the implant recipient coming into the clinic. For more information on Remote Check, please visit <https://www.cochlear.com/us/en/campaign/remote-check-pro>.**

## Conclusion

In conclusion, two years post-FDA approval this paper provides updated guidance based on early experiences to support surgeons and audiologists in the identification, evaluation, counseling, and ongoing management of patients with single-sided deafness or unilateral hearing loss. Through collaboration with leading industry experts, we have identified critical candidacy considerations and emphasized the importance of comprehensive assessments to better understand the quality of life and other impacts on this unique patient population, both pre- and post-implantation.

Understanding patient motivation and establishing realistic expectations are crucial components of this process. These elements ensure that patients are well-prepared for the journey ahead and can achieve the best possible outcomes. Efforts have been made to align evaluation protocols with current MSTB-3 guidelines, and detailed programming guidance has been provided to optimize and streamline the cochlear implant experience.

As we continue to learn and refine these guidelines, there is also an opportunity for Cochlear to develop similar guidance for pediatric patients. This ongoing collaboration aims to ensure optimal outcomes and enhance the overall patient experience, ultimately contributing to the advancement of care for individuals with SSD/UHL.

## Appendix A: Thought leadership meeting participants

Table 1: 2024 Thought Leader Meeting Participants

| Surgeons        | Affiliation   | Audiologists         | Affiliation  |
|-----------------|---|----------------------|--|
| Kevin Brown*    | University of North Carolina<br>School of Medicine<br>Chapel Hill, NC | Allison Beaver       | Rocky Mountain Ear Center<br>Denver, CO                            |
| Daniel Zeitler* | Virginia Mason Franciscan Health<br>Seattle, WA                       | Olivia Boothby       | Virginia Mason Franciscan Health<br>Seattle, WA                    |
| Samantha Anne   | Cleveland Clinic<br>Cleveland, OH                                     | Camille Dunn-Johnson | University of Iowa<br>Iowa City, IA                                |
| Richard Gurgel  | University of Utah<br>Salt Lake City, UT                              | Jill Firszt          | Washington University<br>School of Medicine<br>St. Louis, MO       |
| David Kelsall   | Rocky Mountain Ear Center<br>Englewood, CO                            | Melissa Hurtado      | Mount Sinai<br>New York, NY  |
| Shawn Stevens   | Barrow Neurological Institute<br>Phoenix, AZ                          | Kristen Lewis        | Midwest Ear Institute<br>Kansas City, MO                           |
| Emily Stucken   | University of Michigan<br>Ann Arbor, MI                               | Laurie Mauro         | Children's Hospital of Philadelphia<br>Philadelphia, PA            |
|                 |   | Jessica Novak        | Children's Minnesota<br>Minneapolis, MN                            |
|                 |   | Andrewa Overton      | University of North Carolina School<br>of Medicine Chapel Hill, NC |
|                 |   | William Shapiro      | NYU Langone Health<br>New York, NY                                 |
|                 |   | Molly Smeal          | Cleveland Clinic<br>Cleveland, OH                                  |
|                 |   | Hillary Snapp        | University of Miami Health System<br>Miami, FL                     |
|                 |   | Sarah Sydlowski      | Cleveland Clinic<br>Cleveland, OH                                  |
|                 |   | Viral Tejani         | University Hospitals<br>Cleveland, OH                              |
|                 |   | Gabrielle Watson     | University of Iowa<br>Iowa City, IA                                |
| Psychologist    | Affiliation   |                      |  |
| Ivette Cejas    | University of Miami Miller<br>School of Medicine<br>Miami, FL         |                      |  |

\* Moderator



*Table 2: 2021 Thought Leader Meeting Participants*

| <b>Surgeons</b>      | <b>Affiliation</b>  | <b>Audiologists</b>      | <b>Affiliation</b>   |
|----------------------|---|--------------------------|--|
| Kevin Brown*         | University of North Carolina<br>School of Medicine<br>Chapel Hill, NC | Allison Beiver           | Rocky Mountain Ear Center<br>Denver, CO                      |
| Syed Ahsan           | Kaiser ENT Health<br>Anaheim, CA                                      | Andrea Bucker            | UNC Health<br>Chapel Hill, NC                                |
| Samantha Anne        | Cleveland Clinic<br>Cleveland, OH                                     | Camille Dunn-Johnson     | University of Iowa<br>Iowa City, IA                          |
| Renee Banakis Hartl^ | University of Utah Health<br>Salt Lake City, UT                       | Jill Firszt              | Washington University<br>School of Medicine<br>St. Louis, MO |
| Maura Cosetti        | Mount Sinai<br>New York, NY   | Meredith Holcomb         | University of Miami Health System<br>Miami, FL               |
| Richard Gurgel       | University of Utah<br>Salt Lake City, UT                              | Jourdan Holder           | Vanderbilt University<br>Nashville, TN                       |
| Michael Hoffer^      | University of Miami Health System<br>Miami, FL                        | Laura Schadt             | Baylor College of Medicine<br>Houston, TX                    |
| Jacob Hunter         | UT Southwestern Medical Center<br>Dallas, TX                          | Molly Smeal              | Cleveland Clinic<br>Cleveland, OH                            |
| David Kelsall        | Rocky Mountain Ear Center<br>Denver, CO                               | Johanna Whitson          | UT Southwestern Medical Center<br>Dallas, TX                 |
| Eric Lupo            | Rocky Mountain Ear Center<br>Denver, CO                               |                          |  |
| Brendan O'Connell    | Charlotte Eye, Ear, Nose & Throat<br>Assoc<br>Charlotte, NC           |                          |  |
| Alex Sweeney         | Baylor College of Medicine<br>Otolaryngology<br>Houston, TX           |                          |  |
| Christopher Welch^   | University of Michigan<br>Ann Arbor, MI                               |                          |  |
| Daniel Zeitler       | Virginia Mason Franciscan Health<br>Seattle, WA                       |                          |  |
|                      |   | <b>Hearing Scientist</b> | <b>Affiliation</b>   |
|                      |   | Mario Svirsky            | NYU Langone<br>New York, NY                                  |

\* Moderator    ^ Virtual Attendee

---

## Appendix B: Evidence supporting CI for SSD

Individuals with SSD struggle to process timing and amplitude cues, which are critical for sound localization and speech perception in noisy environments (Gordon et al., 2015). These challenges can lead to listening fatigue, social isolation, and decreased quality of life (Lucas et al, 2018). This may be especially true in complex listening environments such as work, school, or other social situations. In some cases, psychological distress related to fear of losing hearing in the better ear may also impact quality of life (Vannson et al, 2015).

Binaural hearing offers numerous benefits, including improved sound localization, increased loudness perception due to binaural summation, and overall enhanced hearing in both quiet and noisy environments, which can only be achieved with auditory inputs from both ears (Vila & Lieu, 2015). Cochlear implants can help restore some degree of binaural hearing, thereby improving spatial awareness and reducing the cognitive load associated with listening. Early intervention for both adults and children is crucial, as it can prevent the brain from adapting to the loss of input from the deaf ear, thereby preserving binaural auditory processing capabilities. This early intervention ensures that individuals can fully benefit from the advantages of binaural hearing, leading to better auditory outcomes and an improved quality of life.

Tinnitus is also prevalent in this population with some studies reporting up to 86% of individuals with SSD or UHL presenting with some degree of tinnitus (Quaranta et al, 2004). The exact cause of tinnitus is not fully understood, but it is believed that reduced or absent auditory input contributes to changes in neural activity along the auditory pathway (Eggermont & Roberts, 2004). The resulting co-occurrence of hearing loss and tinnitus can significantly impact quality-of-life, making it essential to consider tinnitus in the overall treatment plan and evaluation process. Cochlear implantation has been shown to effectively reduce tinnitus severity while wearing the cochlear implant in a high percentage of cases (50-90%) as reported by Idriss et al, 2022; Assouly et al, 2021; Borges et al, 2021; Levy et al, 2019; while only a small number of cochlear implant cases result in exacerbation or no change in tinnitus perception. This provides strong evidence to support consideration of a CI during the evaluation process.

Although this paper focuses on adults, it is crucial to recognize that children with SSD or UHL also face unique challenges. SSD in children has been correlated to speech and language delays, which can lead to decreased academic performance, behavioral problems, and a reduced quality of life (Jin et al., 2014; Fitzpatrick et al., 2019). Research indicates that SSD can promote an abnormal aural preference, and delayed implantation may hinder the benefits of subsequent implantation (Gordon et al., 2015; Jiwani et al., 2016; Propst et al., 2010; Schmithorst et al., 2014). Addressing these issues early is essential to mitigate their impact and improve outcomes for children with SSD/UHL.

In addition to this, a review of retrospective data was completed by Cochlear and further demonstrates with reasonable assurance the safety and effectiveness of cochlear implantation in individuals with SSD. Full details can be found at [https://www.accessdata.fda.gov/cdrh\\_docs/pdf/P970051S205B.pdf](https://www.accessdata.fda.gov/cdrh_docs/pdf/P970051S205B.pdf).

---

## Appendix C: Assessments for QoL, Tinnitus, and Spatial Hearing

In our efforts to match the right candidate with the optimal solution, quality-of-life and handicap profiles offer clinicians valuable insights for evaluating these factors both before and after implantation. The following tools should be considered and implemented when appropriate.

The Cochlear Implant Quality of Life (CIQOL)- 10 and -35 Profile were derived from the full CIQOL and provides a comprehensive assessment of functional abilities in adults with cochlear implants. Either shortened version may be considered; although the CIQOL-10 does not provide domain specific information but rather a general overview of quality-of-life. Developed at the Medical University of South Carolina, the CIQOL-35 measures outcomes for 6 domains (communication, emotional, entertainment, environment, listening effort, and social). Responses are scored and converted for each domain to an outcome measure score (McRacken & Hand, 2019). These scores can then be compared to normative data to help guide pre- and post-implant counseling. For additional information regarding development and use of the CIQOL-35, please refer to the “Publications” section of MUSC CIQOL website: <https://education.musc.edu/CIQOL>.

The Tinnitus Handicap Inventory (THI), from Newman et al. (1996), was developed as a quick and simple way to evaluate the impact of tinnitus to a patient. It consists of 25 questions and should be considered for use during candidacy evaluation if tinnitus is reported. This assessment can then be compared to post-implant THI scores to measure treatment outcome.

The SSQ-12 is a shortened version of the Speech, Spatial, and Qualities of Hearing Scale (SSQ) consisting of 12 items designed to assess hearing abilities in everyday situations to evaluate hearing abilities that are not reflected on the audiogram. It focuses on three main areas: speech perception, spatial hearing, and qualities of hearing (Noble et al., 2013). For patients with SSD, this tool assists clinicians in evaluating their ability to localize sounds, judge distances, and distinguish between simultaneous sounds.

## References

1. Almugathwi, M., Wearden, A., Green, K., Hill-Feltham, P., & Powell, R. (2020). Online support group users' perceptions and experiences of bone-anchored hearing aids (BAHAs): a qualitative study. *International Journal of Audiology*, 59(11), 850–858.
2. Assouly, K.K.S., Heteren, J.A.A., Stokroos, R.J., Stegeman, I., Smit, A.L. Cochlear implantations for patients with tinnitus. 2021. *Progress in Brain Research*, Volume 260, ISSN 0079-6123, Chapter 2.
3. Borges ALF, Duarte PLES, Almeida RBS, Ledesma ALL, Azevedo YJ, Pereira LV, Bahmad F Jr. Cochlear implant and tinnitus-a meta-analysis. *Braz J Otorhinolaryngol*. 2021 May-Jun;87(3):353-365. doi: 10.1016/j.bjorl.2020.11.006. Epub 2020 Dec 8. PMID: 33342697; PMCID: PMC9422519.
4. Cochlear Limited D1741332 Implant MAP Distribution Statistics. March 2020. Data on file.
5. Dillon M. T., Buss E., Rooth M. A., King E. R., Deres E. J., Buchman C. A., Pillsbury H. C., Brown K. D. (2017b). Effect of cochlear implantation on quality of life in adults with unilateral hearing loss. *Audiol Neurotol*, 22, 259–271.
6. Eggermont JJ, Roberts LE. The neuroscience of tinnitus. *Trends Neurosci* 2004;27:676–82
7. Fitzpatrick, E. M., Gaboury, I., Durieux-Smith, A., Coyle, D., Whittingham, J., & Nassrallah, F. (2019). Auditory and language outcomes in children with unilateral hearing loss. *Hearing Research*, 372, 42–51.
8. Gordon K, Henkin Y, Kral A. Asymmetric hearing during development: the aural preference syndrome and treatment options. *Pediatrics*. 2015;136(1):141-53.
9. Idriss, S.A.; Reynard, P.; Marx, M.; Mainguy, A.; Joly, C.-A.; Ionescu, E.C.; Assouly, K.K.S.; Thai-Van, H. Short- and Long-Term Effect of Cochlear Implantation on Disabling Tinnitus in Single-Sided Deafness Patients: A Systematic Review. *J. Clin. Med*. 2022, 11, 5664. <https://doi.org/10.3390/jcm11195664>
10. Jin SH, Liu C, Sladen DP. The effects of aging on speech perception in noise: comparison between normal-hearing and cochlear-implant listeners. *J Am Acad Audiol*. 2014 Jul-Aug;25(7):656-65. doi: 10.3766/jaaa.25.7.4. PMID: 25365368.
11. Jiwani S, Papsin BC, Gordon KA. Early unilateral cochlear implantation promotes mature cortical asymmetries in adolescents who are deaf. *Hum Brain Mapp*. 2016 Jan;37(1):135-52. doi: 10.1002/hbm.23019. Epub 2015 Oct 12. PMID: 26456629; PMCID: PMC6867517.
12. Kay-Rivest, E, Itrace, AL, Golub, JS, and Svirsky, MA. Prevalence of single-sided deafness in the United States. *Laryngoscope*. (2022) 132:1652–6. doi: 10.1002/lary.29941
13. Kumpik, D. P., & King, A. J. (2019). A review of the effects of unilateral hearing loss on spatial hearing. *Hearing Research*, 372, 17–28.
14. Landsberger DM, Svrakic M, Roland JT, Svirsky M. The Relationship Between Insertion Angles, Default Frequency Allocations, and Spiral Ganglion Place Pitch in Cochlear Implants. *Ear and hearing*. 2015;36(5):e207-e213.
15. Levy DA, Lee JA, Nguyen SA, McRackan TR, Meyer TA, Lambert PR. (2020). Cochlear Implantation for Treatment of Tinnitus in Single-sided Deafness: A Systematic Review and Meta-analysis. *Otol Neurotol*, 41.
16. Lucas, L., Katiri, R., & Kitterick, P. T. (2018). The psychological and social consequences of single-sided deafness in adulthood. *International Journal of Audiology*, 57(1), 21–30.
17. McRackan TR, Hand BN; Cochlear Implant Quality of Life Development Consortium; Vellozo CA, Dubno JR. Cochlear Implant Quality of Life (CIQOL): Development of a Profile Instrument (CIQOL-35 Profile) and a Global Measure (CIQOL-10 Global). *J Speech Lang Hear Res*. 2019 Sep 20;62(9):3554-3563. doi: 10.1044/2019\_JSLHR-H-19-0142. Epub 2019 Sep 4. PMID: 31479616; PMCID: PMC6808347.
18. Nassiri, A. M., Wallerius, K. P., Saoji, A. A., Neff, B. A., Driscoll, C. L. W., Carlson, M. L. (2021a). Impact of duration of deafness on speech perception in single-sided deafness cochlear implantation in adults. *Otol Neurotol*, 43, e45–e49.
19. Newman, C. Jacobson, G. Spitzer, J. Development of the Tinnitus Handicap Inventory. *Arch Otolaryngol Head Neck Surg*. 1996;122(2):143-148.
20. Noble, W., Søgaard Jensen, N., Naylor, G., Bhullar, N. & Akeroyd, M. (2013) A short form of the Speech, Spatial and Qualities of Hearing scale suitable for clinical use: The SSQ12. *Int J Audiol* 52(6), 409–412.
21. Propst EJ, Greinwald JH, Schmithorst V. Neuroanatomic Differences in Children With Unilateral Sensorineural Hearing Loss Detected Using Functional Magnetic Resonance Imaging. *Arch Otolaryngol Head Neck Surg*. 2010;136(1):22–26. doi:10.1001/archoto.2009.208.
22. Quaranta N, Wagstaff S, Baguley DM. Tinnitus and cochlear implantation. *Int J Audiol* 2004;43:245–251.
23. Rahne, T., & Plontke, S. K. (2016). Functional result after cochlear implantation in children and adults with single-sided deafness. *Otol Neurotol*, 37, e332–e340.
24. Schmithorst, Vincent & Plante, Elena & Holland, Scott. (2014). Unilateral deafness in children affects development of multi-modal modulation and default mode networks. *Frontiers in human neuroscience*. 8. 164. 10.3389/fnhum.2014.00164.
25. Vannson, N., James, C., Fraysse, B., Strelnikov, K., Barone, P., Deguine, O., & Marx, M. (2015). Quality of life and auditory performance in adults with asymmetric hearing loss. *Audiology and Neurotology*, 20(Suppl. 1), ~38–43.
26. Vila PM, Lieu JE. Asymmetric and unilateral hearing loss in children. *Cell Tissue Res*. 2015;361(1):271-278.
27. Sagi, E, Azadpour M, Neukam, J, Capach NH, Svirsky MA. (2021) Reducing interaural tonotopic mismatch preserves binaural unmasking in cochlear implant simulations of single-sided deafness. *JASA*. 150(4): 2316-2326. <https://doi.org/10.1121/10.0006446>
28. Tóth, TH, Németh, A, Bakó, P, Révész, P, Gerlinger, I, Szanyi I. (2023). Matching the pitch perception of the cochlear implanted ear with the contralateral ear in patients with sing-sided deafness: a novel approach. *European Archives of Oto-Rhino-Laryngology*. 280:4851-4859.

This material is intended for health professionals. If you are a consumer, please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

©Cochlear Limited 2024. All rights reserved. ACE, Advance OffStylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, 콘트오아, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, SoundBand, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidium are either trademarks or registered trademarks of the Cochlear group of companies.

### Cochlear Americas

10350 Park Meadows Drive, Lone Tree, CO 80124, USA  
Tel: +1 303 790 9010 Support: Tel: +1 800 483 3123

### Cochlear Canada Inc.

2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada  
Tel: +1 800 483 3123 Fax: +1 416 972 5083

[www.cochlear.com/us](http://www.cochlear.com/us)