



Connect and protect

Promoting long term cochlea health

Start strong with Slim Modiolar



Decisions at the time of implantation establish the foundation for long-term outcomes. In particular, electrode design and placement play a critical role in shaping the electrode-neural interface (ENI), which underpins hearing performance.




Two main goals for an electrode


Cochlear's Slim Modiolar electrode is intentionally designed to have the closest connection to the spiral ganglion cells while protecting delicate structures in the cochlea.

Evolution of perimodiolar electrodes

Modern perimodiolar (Slim Modiolar) arrays are not the same as legacy perimodiolar arrays (Contour Advance). Slim Modiolar has **60% less volume** than Contour Advance.



Slim Modiolar electrode




Contour Advance® electrode

Connect to the hearing nerve



- 1 Spiral ganglion cells:
Electrode-neuron interface
site for electrical stimulation¹⁻⁷
- 2 Scala vestibuli
- 3 Scala media
- 4 Spiral ligament
- 5 Basilar membrane
- 6 Scala tympani
- 7 **Perimodiolar electrode**
- 8 **Lateral wall electrode**



Focused stimulation and electrode discrimination

Cochlear is the only company to offer **22 electrode contacts**—supporting detailed frequency representation and perimodiolar placement designed to bring electrodes closer to the modiolus.⁷

Perimodiolar placement is associated with significantly improved electrode discrimination success (**77% versus 47%**) and reduced channel interaction.⁸

Superior speech perception outcomes

Clinical outcomes support electrode design as a significant contributor to post-operative speech perception. In a cross-manufacturer study, perimodiolar electrodes demonstrate an average **CNC word score advantage** of approximately **16.9% at 6 months** and **19.9% at 24 months** compared with lateral wall arrays.⁹

When different electrode array types are implanted in the same patient, performance differences can be attributed to the electrode design. In an inpatient comparison, **perimodiolar electrodes averaged 66% CNC word scores compared to a 42% average with lateral wall electrodes.**¹⁰

When controlling for cochlea health, wrapping factor, angular insertion depth, and programming strategy, perimodiolar electrodes with **tighter modiolar proximity** outperform lateral wall electrodes at 6 months (**53.3% versus 45.1%**), emphasizing the importance of these factors in optimizing cochlear implant outcomes.¹¹

Pediatric clinical evidence

Pediatric data indicate that closer modiolar proximity is associated with **lower impedances** and **reduced neural thresholds, reduced energy consumption**, and **more consistent early auditory access**. These factors translated to significantly **improved speech and language** outcomes in children implanted at three years of age or younger.¹²

22 electrode contacts and perimodiolar positioning enhance spectral resolution and maximize speech performance⁷



22 electrode contacts




16 electrode contacts



12 electrode contacts

CNC word scores¹⁰




Protect the cochlea



Cochlear anatomy and insertion trauma

Anatomical studies show that the scala tympani narrows asymmetrically with increasing insertion depth, with lateral wall height decreasing more rapidly than modiolar height beyond approximately 360 degrees. This anatomic change **increases the risk of basilar membrane contact and trauma with deep lateral wall insertions**.^{13,14}




Biological response and fibrosis

Post-implant foreign body response results in new tissue and fibrotic growth regardless of electrode design. Lateral wall electrodes sit closer to the basilar membrane and vascular supply, whereas **perimodiolar electrodes are positioned farther from these vulnerable structures**, potentially reducing long-term impact.²⁰⁻²¹



Perimodiolar electrode

By avoiding contact with the lateral wall of the cochlea, you could potentially **avoid injuring the rich blood supply** which can lead to intracochlear bleeding and inflammation.²⁰




Lateral wall electrode

Even if contact between a straight electrode and basilar membrane can be avoided, the fibrous tissue formed around the electrode array may **impede basilar membrane function** after surgery.²¹

Cochlear duct length ≠ spiral ganglion coverage

Spiral ganglion neurons terminate before the full length of the cochlear duct (average 695 degrees)¹⁵, meaning **deeper insertions may not provide additional neural benefit**.



Slim Modiolar electrodes have **lower translocation rates** compared to legacy perimodiolar electrodes and translocate less than lateral wall electrodes.^{19,*}

Shaul, 2020¹⁶

0.0% vs. 17%

CI532 vs. CI512

Beck, 2023¹⁷

2.7%


CI532 or CI632

Liebscher, 2021¹⁸

5.1% vs. 32.3%

CI532 vs. CI512

Translocation rate*



Hearing preservation and long-term stability

A meta-analysis published by Geng et al. concluded that slim modiolar electrodes (SME) and slim straight electrodes (SSE) have similar hearing preservation abilities in the short term. However, a prominent decline in hearing preservation over time was observed in the SSE compared with SME. This indicates SME may offer better durability in HP over time.²²

Summary

Global insights on the electrode-neural interface

An international group of CI surgeons each with more than 30,000 patients in their care provide a collective expert opinion narrative summarizing the identified influencing factors upon long-term CI outcomes.²³



Connect

“The proximity of the electrode contacts to the auditory nerve seems to be beneficial in terms of resolution and ultimately CI performance.”




Protect

“Perimodiolar positioning causes a much lower risk of disrupting lateral wall structures, thus better preserving long-term cochlear health.”

Clinical takeaway

The Slim Modiolar electrode is engineered to achieve consistent perimodiolar positioning while minimizing intracochlear trauma, thereby optimizing the relationship between electrode contacts and spiral ganglion neurons. Current literature supports the position that electrodes achieving closer modiolar proximity and avoiding delicate cochlear structures establish more favorable conditions for neural stimulation and support both structural preservation and functional outcomes. These findings underscore the importance of electrode choice at the time of implantation.

To learn more about our Slim Modiolar electrode, visit [**cochlear.us/connectandprotect**](https://cochlear.us/connectandprotect)



References

- Goldwyn JH, Bierer SM, Bierer JA. Modeling the electrode-neuron interface of cochlear implants: effects of neural survival, electrode placement, and the partial tripolar configuration. *Hear Res*. 2010 Sep 1;268(1-2):93-104. doi: 10.1016/j.heares.2010.05.005. Epub 2010 May 24.
- Saunders, E., Cohen, L., Aschendorff, A., Shapiro, W., Knight, M., Stecker, M., Richter, B., Waltzman, S., Tykocinski, M., Roland, T., Laszig, R., & Cowan, R. (2002). Threshold, comfortable level and impedance changes as a function of electrode-modiolar distance. *Ear and hearing*, 23(1 Suppl), 28S-40S. https://doi.org/10.1097/00003446-200202001-00004.
- Cohen, L. T., Richardson, L. M., Saunders, E., & Cowan, R. S. (2003). Spatial spread of neural excitation in cochlear implant recipients: comparison of improved ECAP method and psychophysical forward masking. *Hearing research*, 179(1-2), 72-87. https://doi.org/10.1016/s0378-5955(03)00096-0.
- Hughes, M. L., & Abbas, P. J. (2006). Electrophysiologic channel interaction, electrode pitch ranking, and behavioral threshold in straight versus perimodiolar cochlear implant electrode arrays. *The Journal of the Acoustical Society of America*, 119(3), 1538-1547. https://doi.org/10.1121/1.2164969.
- Basta, D., Todt, I., & Ernst, A. (2010). Audiological outcome of the pull-back technique in cochlear implant recipients. *The Laryngoscope*, 120(7), 1391-1396. https://doi.org/10.1002/lary.20942.
- DeVries, L., Scheperle, R., & Bierer, J. A. (2016). Assessing the Electrode-Neuron Interface with the Electrically Evoked Compound Action Potential, Electrode Position, and Behavioral Thresholds. *Journal of the Association for Research in Otolaryngology : JARO*, 17(3), 237-252. https://doi.org/10.1007/s10162-016-0557-9.
- Berg, K. A., Noble, J. H., Dawant, B. M., Dwyer, R. T., Labadie, R. F., & Gifford, R. H. (2019). Speech recognition as a function of the number of channels in perimodiolar electrode recipients. *The Journal of the Acoustical Society of America*, 145(3), 1556. https://doi.org/10.1121/1.5092350.
- Ramos de Miguel, Á., Argudo, A. A., Borkoski Barreiro, S. A., Falcón González, J. C., & Ramos Macías, A. (2018). Imaging evaluation of electrode placement and effect on electrode discrimination on different cochlear implant electrode arrays. *European archives of oto-rhino-laryngology : official journal of the European Federation of Oto-Rhino-Laryngological Societies (EUFOS) : affiliated with the German Society for Oto-Rhino-Laryngology - Head and Neck Surgery*, 275(6), 1385-1394. https://doi.org/10.1007/s00405-018-4943-2.
- Pennington-FitzGerald, W., Waring, N. A., Hamiter, M., Kuhlmeier, M., & Kim, A. H. (2025). Impact of Cochlear Implant Electrode Array Design on Post-Op Speech Perception. *Otolaryngology--head and neck surgery : official journal of American Academy of Otolaryngology-Head and Neck Surgery*, 172(3), 960-966. https://doi.org/10.1002/ohn.995.
- Ewer, N., Espahbodi, M., Durham, A. R., Johnson, K., Gurgel, R. K., & Patel, N. S. (2026). A Unique Perspective: Inpatient Comparison of Perimodiolar and Lateral Wall Cochlear Implant Electrodes. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 10.1097/MAO.0000000000004860. Advance online publication. https://doi.org/10.1097/MAO.0000000000004860.
- Walia, A., Shew, M. A., Smetak, M., Ortmann, A., Durakovic, N., Jiramongkolchai, P., Leffer, S., Herzog, J. A., & Buchman, C. A. (2026). Perimodiolar Electrode Locations Outperform Lateral Wall Arrays When Controlling for Cochlear Health and Speech Processing Strategy. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 47(2), 257-265. https://doi.org/10.1097/MAO.0000000000004699.
- Tuset, M. P., Bois, E., Noel-Petroff, N., Teissier, N., Van Den Abbeele, T., & Benoit, C. (2026). Comparison of Slim Perimodiolar and Slim Straight Electrode Arrays in Pediatric Cochlear Implantation: A Radiologic, Electrophysiological, and Audiological Study. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 47(3), 447-455. https://doi.org/10.1097/MAO.0000000000004784.
- Avci, E., Nauwelaers, T., Lenarz, T., Hamacher, V., & Kral, A. (2014). Variations in microanatomy of the human cochlea. *The Journal of comparative neurology*, 522(14), 3245-3261. https://doi.org/10.1002/cne.23594.
- Rebscher, S. J., Hetherington, A., Bonham, B., Wardrop, P., Whinney, D., & Leake, P. A. (2008). Considerations for design of future cochlear implant electrode arrays: electrode array stiffness, size, and depth of insertion. *Journal of rehabilitation research and development*, 45(5), 731-747. https://doi.org/10.1682/jrrd.2007.08.0119.
- Danielian, A., Ishiyama, G., Lopez, I. A., & Ishiyama, A. (2020). Morphometric linear and angular measurements of the human cochlea in implant patients using 3-dimensional reconstruction. *Hearing research*, 386, 107874. https://doi.org/10.1016/j.heares.2019.107874.
- Shaul, C., Weder, S., Tari, S., Gerard, J. M., O'Leary, S. J., & Briggs, R. J. (2020). Slim, Modiolar Cochlear Implant Electrode: Melbourne Experience and Comparison With the Contour Perimodiolar Electrode. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 41(5), 639-643. https://doi.org/10.1097/MAO.0000000000002617.
- Beck, R., Aschendorff, A., Arndt, S., Hildenbrand, T., Rauch, A. K., & Ketterer, M. C. (2024). Evaluation of insertion quality of a slim perimodiolar electrode array. *European archives of oto-rhino-laryngology : official journal of the European Federation of Oto-Rhino-Laryngological Societies (EUFOS) : affiliated with the German Society for Oto-Rhino-Laryngology - Head and Neck Surgery*, 281(3), 1215-1220. https://doi.org/10.1007/s00405-023-08212-5.
- Liebscher, T., Mewes, A., Hoppe, U., Hornung, J., Brademann, G., & Hey, M. (2021). Electrode Translocations in Perimodiolar Cochlear Implant Electrodes: Audiological and Electrophysiological Outcome. *Zeitschrift für medizinische Physik*, 31(3), 265-275. https://doi.org/10.1016/j.zemedi.2020.05.004.
- Dhanasingh, A., & Jolly, C. (2019). Review on cochlear implant electrode array tip fold-over and scalar deviation. *Journal of otology*, 14(3), 94-100. https://doi.org/10.1016/j.joto.2019.01.002.
- Bester, C, et al. Four-point impedance as a biomarker for bleeding during cochlear implantation. *Scientific Reports*. 2020;10(1). Razmovski, T., Bester, C., Collins, A., Tan, E., & O'Leary, S. J. (2022). Four-Point Impedance Changes After Cochlear Implantation for Lateral Wall and Perimodiolar Implants. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 43(10), e1107-e1114. https://doi.org/10.1097/MAO.0000000000003732.
- Verberne J, et al. The Effect of Scala Tympani Morphology on Basilar Membrane Contact With a Straight Electrode Array: A Human Temporal Bone Study. *Otol Neurotol*. 2017 Jan;38(1):47-53.
- Geng, S., Wijewickrema, S., Copson, B., Gerard, J. M., & O'Leary, S. (2025). Hearing Preservation of Slim Modiolar and Slim Straight Electrodes: A Systematic Review and Meta-Analysis. *Otology & neurotology : official publication of the American Otological Society, American Neurotology Society [and] European Academy of Otolology and Neurotology*, 46(7), 733-742. https://doi.org/10.1097/MAO.0000000000004540.
- Ramos Macías, A., Ramos-Macías, Á., Briggs, R., Choi, B. Y., Friedmann, D., Ishiyama, A., Lenarz, T., Mylanus, E., O Apos Leary, S., Roland, J. T., Jr, & Zarowski, A. (2025). The Importance of the Electrode-Neural Interface in Supporting Long-Term Outcomes in Cochlear Implantation: Expert Opinion. *Audiology & neuro-otology*, 30(5), 459-470. https://doi.org/10.1159/000546003.

* Cross study comparisons are limited by differences in study design, populations, and endpoints.

This material is intended for health professionals. If you are a consumer, please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

The content of this guideline is intended as a guide for information purposes only and does not replace or remove clinical judgment or the professional care and duty necessary for each specific recipient case. Clinical care carried out in accordance with this guideline should be provided within the context of locally available resources, expertise, and standards of care and practice.

©Cochlear Limited 2026. All rights reserved. ACE, Advance OffStylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, Nexa, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, SoundBand, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidum are either trademarks or registered trademarks of the Cochlear group of companies.

